

ORDERING GUIDE

ELECTROSURGERY FORCEPS

FORCEPS

STANDARD BIPOLAR FORCEPS



Ref.

ESU Plug: 2 pin standard . 3 m cable. Single Use.

FR-700-1105C	McPherson Bipolar Forceps, Straight, 0.5mm Tip, 10.8cm Length.	10u
FR-720-1210C	Adson Bipolar Forceps, 1.0mm Tip, 12.1cm Length.	10u
FR-730-1505C	Cushing Bipolar Forceps, 0.5mm Tip, 15cm Length.	10u
FR-730-1510C	Straight Bipolar Forceps, 1.0mm Tip, 15cm Length.	10u
FR-730-1810C	Straight Bipolar Forceps, 1.0mm Tip, 18cm Length.	10u
FR-730-2010C	Straight Bipolar Forceps, 1.0mm Tip, 20cm Length.	10u
FR-731-2010C	Straight Bipolar Forceps, 1.0mm Angled Tip, 20cm Length.	10u
FR-750-1105C	Jeweler n.5 Bipolar Forceps, Straight, 0.5mm Tip, 11.4cm Length.	10u
FR-760-1610C	Jensen Bayonet Bipolar Forceps, 1.0mm Tip, 16.5cm Length.	10u
FR-760-2010C	Jensen Bayonet Bipolar Forceps, 1.0mm Tip, 20cm Length.	10u
FR-760-2210C	Jensen Bayonet Bipolar Forceps, 1.0mm Tip, 22cm Length.	10u
FR-760-2410C	Jensen Bayonet Bipolar Forceps, 1.0mm Tip, 24cm Length.	10u

FORCEPS

STANDARD MONOPOLAR FORCEPS



Ref.

ESU Plug: 3 pin standard*. 3 m cable. Single Use.

FR-380-1810C	Potts-Smith Monopolar forceps length 7" (17.8cm) 1.0mm tip serrated.	10u
FR-380-1815C	Potts-Smith Monopolar forceps length 7" (17.8cm) 1.5mm tip serrated.	10u
FR-380-2010C	Potts-Smith Monopolar forceps length 8" (20.3cm) 1.0mm tip serrated.	10u
FR-380-2015C	Potts-Smith Monopolar forceps length 8" (20.3cm) 1.5mm tip serrated.	10u
FR-380-2410C	Potts-Smith Monopolar forceps length 9.5" (24cm) 1.0mm tip serrated.	10u
FR-380-2415C	Potts-Smith Monopolar forceps length 9.5" (24cm) 1.5mm tip serrated.	10u
FR-390-2015C	Jansen Bayonet Monopolar forceps, length 8" (20.3cm) 1.5mm tip serrated.	10u