

ORDERING GUIDE

ELECTROSURGERY PENCIL

PENCIL

STANDARD PENCIL WITH BLADE ELECTRODE



Ref.	NONSTICK™ BLADE ELECTRODE. Monopolar. ESU Plug: 3 pin Standard*. 3 m cable. Single Use.	
PMS252	Push buttons Pencil with 70 mm NonStick™ blade electrode.	25u
PMS254	Push buttons Pencil with 70 mm NonStick™ blade electrode, Shielded tip.	25u
PMS252NBH	Push buttons Pencil with 70 mm NonStick™ blade electrode, with Holster.	25u
PMS254NBH	Push buttons Pencil with 70 mm NonStick™ blade electrode, Shielded tip, with Holster.	25u
PMS252NRH	Rocker switch Pencil with 70 mm NonStick™ blade electrode, with Holster.	25u
PMS254NRH	Rocker switch Pencil with 70 mm NonStick™ blade electrode, Shielded tip, with Holster.	25u
PMS248H	Footswitch activated Pencil with 70 mm NonStick™ blade electrode, with Holster.	25u
Ref.	STAINLESS STEEL BLADE ELECTRODE. Monopolar. ESU Plug: 3 pin Standard*. 3 m cable. Single Use.	
PMS250H	Push buttons Pencil with 70 mm Stainless Steel blade electrode, with Holster.	25u
PMS249H	Footswitch activated Pencil with 70 mm Stainless Steel blade electrode, with Holster	25u