

## ORDERING GUIDE

## MACROSEAL

## MACROSEAL

## MONOPOLAR FORCEPS



## Ref.

**HAND-ACTIVATED. ESU Plug: 3 pin standard\*. 3m tube. 4m cable. Single Use.**

MS8805	Hand-activated monopolar smoke evacuation forceps, 0.5 mm NonStick™ tips.	20u
MS8808	Hand-activated monopolar smoke evacuation forceps, 0.75 mm NonStick™ tips.	20u
MS8812	Hand-activated monopolar smoke evacuation forceps, 1.2 mm NonStick™ tips.	20u
MS8812SS	Hand-activated monopolar smoke evacuation forceps, 1.2 mm SS tips.	20u

## Ref.

**FOOTSWITCH-ACTIVATED. ESU Plug: 3 pin standard\*. 3m tube. 4m cable. Single Use.**

MF8805	Footswitch-activated monopolar smoke evacuation forceps, 0.5 mm NonStick™ tips.	20u
MF8808	Footswitch-activated monopolar smoke evacuation forceps, 0.75 mm NonStick™ tips.	20u
MF8812	Footswitch-activated monopolar smoke evacuation forceps, 1.2 mm NonStick™ tips.	20u
MF8812SS	Footswitch-activated monopolar smoke evacuation forceps, 1.2 mm SS tips.	20u

## Ref.

**HAND-ACTIVATED. ESU Plug: 3 pin standard\*. 3m cable. Single Use.**

MS8605	Hand-activated monopolar forceps, 0.5 mm NonStick™ tips (WITHOUT smoke evacuation).	25u
MS8608	Hand-activated monopolar forceps, 0.75 mm NonStick™ tips (WITHOUT smoke evacuation).	25u
MS8612	Hand-activated monopolar forceps, 1.2 mm NonStick™ tips (WITHOUT smoke evacuation).	25u
MS8612SS	Hand-activated monopolar forceps, 1.2 mm SS tips (WITHOUT smoke evacuation).	25u