

ELECTROSURGERY

PENCILS FORCEPS & CABLES



PRIMA

COMPACT PORTFOLIO CATALOGUE

SUPERVAC

TELESCOPIC SMOKE EVACUATION PENCIL

Exclusive sliding buttons allow telescopic extension from 90-150mm.

- Powerful evacuation performance – max flow rate of 85L/min.
- Provides greater visualisation of the surgical site, especially in cavities.
- 360° Swivel connection and lightweight construction enable freedom of movement, reducing wrist fatigue.

Available in 5 configurations:



ACTISEAL

HAND-ACTIVATED BIPOLAR SMOKE EVACUATION FORCEPS

The world's first and only forceps to offer integrated hand-activation and smoke evacuation.

Actiseal™ is ideally suited to fine, low-powered superficial procedures.

PRECISION AT YOUR FINGERTIPS

Actiseal "two-stage" buttons allows to easily:

- 1) Manipulate tissue without activation
- 2) Greater pressure activates the current

NO FOOTSWITCH REQUIRED

MACROSEAL

HAND-ACTIVATED MONOPOLAR SMOKE EVACUATION FORCEPS

The world's first and only forceps to offer hand-activation and integrated smoke evacuation, enabling the precise manipulation, transection and coagulation of target tissues at the touch of a button.




The "buzz me" technique, responsible for the burnt fingers of surgeons everywhere, can be completely eliminated.

NO FOOTSWITCH REQUIRED



Squeeze the buttons to grip tissue

Squeeze harder until "click" to activate power



ClearAir²

SMOKE EVACUATION SYSTEM

SAFER


ClearAir²™ utilises powerful ULPA filtration to remove 99.999% of harmful particulates whilst monitoring filter status.

SMARTER

ClearAir²™ turns on/off automatically without the need for additional cables or accessories.

QUIETER

Acoustically engineered to absorb vibrations and dampen noise, providing for greater concentration and comfort.



ClearAir² Unit



CLEAR Filter

Daily ULPA filter

GREEN Filter

Weekly ULPA filter

RETRACT


MICRODISSECTION PENCIL

Specifically designed for microsurgery, Retract™ contains a ultra-sharp tungsten microdissection needle which can be extended and retracted.

THE SURGEON'S BALL-POINT PEN

It's refined ergonomics allow for an incredible level of precision and provides for absolute control in very delicate areas.

PRECISION. SAFETY. CONTROL.



ELECTROSURGERY PENCIL

WITH NONSTICK™ BLADE

Supplied with a choice of 3 different control schemes and electrode configurations:


- A** 70 mm Nonstick™ Blade
- B** 70 mm Nonstick™ Blade, Shielded Tip
- C** 70mm Stainless Steel Blade

Nonstick™ electrodes reduce eschar adhesion and improve cutting & coagulation performance.

Available with safety holster



No scratch pad required



ORDERING GUIDE

SMOKE EVACUATION PORTFOLIO

REF.	PRODUCT DESCRIPTION	MONOPOLAR	BIPOLAR	STD 3 PIN	* 3 PIN	** 2 PIN	4 m	10/22 mm	3 m	QTY
SUPERVAC	TELESCOPIC PENCIL									
PMS6300N	With 70 mm NONSTICK Blade Electrode	X	-	X	-	-	X	X	X	20
PMS6305N	With 70 mm NONSTICK Blade Electrode, Shielded Tip	X	-	X	-	-	X	X	X	20
PMS6303N	With 70 mm Ceramic Blade Electrode	X	-	X	-	-	X	X	X	20
PMS6301N	With 70 mm Tungsten Needle Electrode	X	-	X	-	-	X	X	X	20
PMS6302N	Without Electrode	X	-	X	-	-	X	X	X	20
ACTISEAL	HAND-ACTIVATED BIPOLAR FORCEPS									
AS8805	0.5 mm Tips	-	X	-	X	-	X	X	X	20
AS8808	0.75 mm Tips	-	X	-	X	-	X	X	X	20
AS8812	1.2 mm Tips	-	X	-	X	-	X	X	X	20
ACTISEAL	FOOTSWITCH-ACTIVATED BIPOLAR FORCEPS									
BF8805	0.5 mm Tips	-	X	-	-	X	X	X	X	20
BF8808	0.75 mm Tips	-	X	-	-	X	X	X	X	20
BF8812	1.2 mm Tips	-	X	-	-	X	X	X	X	20
MACROSEAL	HAND-ACTIVATED MONOPOLAR FORCEPS									
MS8805	0.5 mm NONSTICK Tips	X	-	X	-	-	X	X	X	20
MS8808	0.75 mm NONSTICK Tips	X	-	X	-	-	X	X	X	20
MS8812	1.2 mm NONSTICK Tips	X	-	X	-	-	X	X	X	20
MS8812SS	1.2 mm Stainless Steel Tips	X	-	X	-	-	X	X	X	20
MACROSEAL	FOOTSWITCH-ACTIVATED MONOPOLAR FORCEPS									
MF8805	0.5 mm NONSTICK Tips	X	-	X	-	-	X	X	X	20
MF8808	0.75 mm NONSTICK Tips	X	-	X	-	-	X	X	X	20
MF8812	1.2 mm NONSTICK Tips	X	-	X	-	-	X	X	X	20
MF8812SS	1.2 mm Stainless Steel Tips	X	-	X	-	-	X	X	X	20
CLEARAIR²	SURGICAL SMOKE EVACUATION SYSTEM									
PMS900SE	ClearAir ² Unit	X	X	-	-	-	-	-	-	1
CFU-001	CLEAR Replacement ULPA Filter For ClearAir ² Unit	X	X	-	-	-	-	-	-	10
GFU-001	GREEN Replacement ULPA Filter For ClearAir ² Unit	X	X	-	-	-	-	-	-	4

NON-SMOKE EVACUATION PORTFOLIO

REF.	PRODUCT DESCRIPTION	MONOPOLAR	BIPOLAR	STD 3 PIN	4 mm BANANA PLUG	* 3 PIN	3 m	QTY
RETRACT	MICRODISSECTION PENCIL							
PMS8001	With Adjustable (7mm/14mm) / Retractable Needle	X	-	X	-	-	X	10
PENCIL	With 70mm NONSTICK™ BLADE ELECTRODE							
PMS252	Push Buttons	X	-	X	-	-	X	25
PMS254	Push Buttons, Shielded Tip Blade Electrode	X	-	X	-	-	X	25
PMS252NBH	Push Buttons + Holster	X	-	X	-	-	X	25
PMS254NBH	Push Buttons, Shielded Tip Blade Electrode + Holster	X	-	X	-	-	X	25
PMS252NRH	Rocker Switch + Holster	X	-	X	-	-	X	25
PMS254NRH	Rocker Switch, Shielded Tip Blade Electrode + Holster	X	-	X	-	-	X	25
PMS248H	Footswitch Activated + Holster	X	-	-	X	-	X	25
	With 70mm STAINLESS STEEL BLADE ELECTRODE							
PMS249H	Footswitch Activated + Holster	X	-	-	X	-	X	25
PMS250H	Push Buttons + Holster	X	-	X	-	-	X	25
ACTISEAL	HAND-ACTIVATED BIPOLAR FORCEPS							
AS8605	0.5 mm Stainless Steel Tips	-	X	-	-	X	X	25
AS8608	0.75 mm Stainless Steel Tips	-	X	-	-	X	X	25
AS8612	1.2 mm Stainless Steel Tips	-	X	-	-	X	X	25
MACROSEAL	HAND-ACTIVATED MONOPOLAR FORCEPS							
MS8605	0.5 mm NONSTICK Tips	X	-	X	-	-	X	25
MS8608	0.75 mm NONSTICK Tips	X	-	X	-	-	X	25
MS8612	1.2 mm NONSTICK Tips	X	-	X	-	-	X	25
MS8612SS	1.2 mm Stainless Steel Tips	X	-	X	-	-	X	25

* 3 pin plug for use with Valleylab, Megadyne, Conmed and Bovie electrosurgery generators







** 2 pin plug for use with Erbe, Bowa and KLS Martin electrosurgery generators




UNIVERSAL BIPOLAR CABLE

2in1 Fits US and European bipolar instruments.
Reusable up to 100 autoclave cycles.

ORDERING INFORMATION

REF.	PRODUCT DESCRIPTION	PLUG		
P-30	Reusable cable for bipolar instruments. Plug for Erbe, Siemens Sutter BM50.		4 m	1 u.
Q-30	Reusable cable for bipolar instruments. Plug with 2 mm pins for Martin, Aesculap, Integra, Wolf, Berchtold, Sutter BM780.		4 m	1 u.
S-30	Reusable cable for bipolar instruments. Plug with twin 4 mm pin, Erbe, Eschmann, Alsa MCD series.		4 m	1 u.
V-30	Reusable cable for bipolar instruments. Block connector (28mm), 4 mm pins for Bowa, EMC, Erbe, Ethicon, Olympus, Soering, Conmed 5000.		4 m	1 u.





UNIVERSAL MONOPOLAR CABLE

2in1 4 / 8 mm Male plug & plug Female plug for 4-5 mm pin.
Reusable up to 100 autoclave cycles.



ORDERING INFORMATION

REF.	PRODUCT DESCRIPTION		
PMS8700	Monopolar Universal, 2 in 1 generator connector (8 mm / 4 mm) to Universal Female plug.	4 m	1 u.

This document is intended solely for the use of healthcare professionals according to his/her own professional clinical judgment. The information presented is only intended to demonstrate Prima's products, not to recommend a particular technique or procedure.

Products may not be available in all markets due to regulatory or commercial reasons.

No part of this document can be photocopied, reproduced or translated into any other languages without written consent.

Prima logo and products, are registered trademarks of Prima Medical Ltd. All other trademarks are trademarks of their respective owners or holders. The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Prima's trademark or other intellectual property rights concerning that name or logo. All rights reserved.

PRIMA
ELECTROSURGERY INNOVATION

Prima Medical Limited
The Counting House,
Bond's Mill Estate,
Stonehouse, GL10 3RF, UK

P: +44 (0) 1453 707671 (ext2)
M: +44 (0) 7539 662635

info@prima-medical.com
www.prima-medical.com